DTR.M.005 FRAZIER SUCTION HANDLES

Instructions for Use

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DESCRIPTION

Frazier Suction Devices can be used in surgical procedures where fine, delicate suctioning is necessary to evacuate blood or other fluids from confined spaces or from the surgical site. Some Frazier Suction Devices are controlled by a small opening on the handle. Other variations include fenestration, cross fenestration, straight and bent to fit the user needs.

MATERIALS OF CONSTRUCTION

Metal: FSH2006, FSH2006-B, FSF2006, FSF2046, FSN2426, FSH2008, FSH2008-B, FSF2008, FSF2048, FSH2009, FSH2009-B, FSF2009, FSF2049, FSN2429, FSH2012, FSF2012, FSF2412, FSN24212 and FSH2015

Product Component	Details
Metal Tube	316L Stainless Steel
Frazier Head/Handle	303/304 Stainless Steel
Stylet (not included with all Frazier products)	302 Stainless Steel

Metal and Plastic: 74-2017 and 74-2050

Product Component	Details
Metal Tube	304 Stainless Steel
Frazier Head/Handle	ABS Plastic Hub

INTENDED USE/INDICATIONS FOR USE

The Frazier Suction Handles are indicated for the use of aspiration of blood and other bodily fluids from the surgical or diagnostic operation in ENT procedures.

INTENDED USER

The intended user of the Frazier Suction Handle consists of qualified healthcare professionals and medical experts who have been trained in suction procedures, and in the use of aspirators in a hospital / clinical environment.

CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

There are no absolute contraindications but there are risks associated with suctioning and this should be weighed against patient specific needs.

PATIENT POPULATION

No age restriction. As the devices were developed to support suction activities during a surgical procedure, the target population for this device is anyone, of any age, in the human population, that is identified as needing suction during an ENT or General Surgery procedure.

WARNINGS AND PRECAUTIONS

- Device is single use only.
- Device is provided sterile, in a single sterile barrier packaging. Do not sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- Do not use if the packaging is damaged. Packaging should be carefully inspected for punctures, tears, security of the seal, and any other evidence of the sealed pouch having been compromised prior to the placement of the contents in the sterile field.
- Device should be used for suction only – it is not suitable for irrigation uses.
- Additional care should be taken if patients are suffering from bradycardia, hypoxia or other similar complaints / concerns.
- There is a risk of epistaxis if patient moves or struggles while the device is in use.
- Excessive suction can cause tissue damage.
- Suctioning is not entirely without risk, as with all medical interventions, healthcare workers must use clinical judgement to determine whether the benefit outweighs the risks.
- Do not use in the presence of magnetic fields such as those which occur during magnetic resonance imaging. The device has not been tested for magnetically induced torque, radiofrequency induced heating or image artefacts.
- Instruments made from Stainless Steel may contain Nickel. Although the contact between instruments and patient is transient, it may cause

an allergic reaction in patients with sensitivity to this material. Users should risk assess the use of the instruments in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.

SHELF LIFE AND STERILITY

- Instruments are supplied in a sterile condition and are packaged individually. The instruments are sterilised by exposure to ethylene oxide indicated by the  symbol on the label. They remain sterile as long as the package integrity has not been violated.
- Packaging must be inspected before use. Do not use any component from an opened or damaged package.
- The expiry date is printed on the label. Do not use instruments after the expiry date.

STORAGE AND HANDLING

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight, and temperature/humidity extremes. Store in a cool, dry place.
- Care must be exercised in the handling of cases or individual instruments to prevent damage to the sterile barrier.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.


POSSIBLE ADVERSE EFFECTS


- Complications such as discomfort, bronchoconstriction, infection, injury to the tracheal mucosa, hemorrhage, atelectasis and hemodynamic changes could occur as a result of the procedure.
- The most frequently reported complications are reduced oxygenation and pulmonary de-recruitment. Hyper-oxygenating the patient before suctioning and allowing rest between suctioning attempts can reduce the chances of hypoxia.
- Injury to the mucosal membranes and bleeding can also occur.

- Using the appropriate vacuum level, a water-soluble lubricant and avoiding forcing the device can reduce the risk of tissue injury.

OPERATING INSTRUCTIONS

1. Aspirating procedures can begin after it has been explained to the patient and supplies have been collected. The use of PPE is highly recommended and should be in line with hospital/clinic processes. Sterile technique should always be used in order to reduce the risk of nosocomial infections.
2. Frazier Suction Handles are intended to be used with silicone tubing and an aspirator. A stylet may also be used to remove blockages of the Suction Handle during use. The Stylet may be provided with the product, in the sterile packaging, and should be kept in a suitable place until it is required, minimizing any contamination or de-sterilisation.
3. The metal Frazier Suction Handles may be bent to fit procedure requirements. This should not exceed an angle of 40 degrees (shown below) and the recommended distances below should be followed. Please note that manipulation of the device may adversely affect the suction of the device.

Bending points at the minimum distance.

DEVICE DISPOSAL

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste. Handle and dispose this product in accordance with accepted medical practice and with applicable local, state and national laws and regulations.

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC