DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY

CERVICAL ROTATING BIOPSY PUNCH

Instructions for Use

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REF CBR1041


DESCRIPTION

The Cervical Rotating Biopsy Punch has been designed specifically for use in Gynaecological procedures to obtain a cervical biopsy. Equipped with a pair of jaws at the end of a metal rod with a handle and trigger. The long metal construction ensures better access to the cervix for the jaws, whilst the handle is operated outside.

When the trigger is pulled the sharp jaw closes, this action removes a segment of tissue for the biopsy. 360° rotation enables enhanced positioning and patient interaction.

MATERIALS OF CONSTRUCTION

Product Component	Details
Top Jaw	17-4 Precipitation Hardened Stainless Steel
Bottom Jaw	316 Stainless Steel
Tube	304 Stainless Steel
Handle	Polycarbonate (PC)
Trigger	Polycarbonate (PC)
Rotation Controller	Polycarbonate (PC)

INTENDED USE/INDICATIONS FOR USE

The Cervical Rotating Biopsy Punch is indicated for use in Gynaecological procedures to obtain a cervical biopsy and allows a biopsy to be taken from a cervix.

INTENDED USER

The Cervical Rotating Biopsy Punch is intended for use by Healthcare Practitioners trained in Biopsy sampling procedures. It must only be used by suitably trained and qualified staff. The product is intended solely for use in the medical sector and within clinical environment. The user shall be familiar with the instrument before using it.

INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

The Cervical Rotating Biopsy Punch is indicated for patients requiring cervical biopsy.

The Cervical Rotating Biopsy Punch is to be used as per the clinician’s decision based on cervical finding or suspicions of related to medical indications such as but not limited to:

- Abnormal PAP smear (HPV)
- Lesions of the cervix
- Atypical squamous cells of undetermined significance (ASCUS)
- Low-grade squamous intraepithelial lesions of the cervix (LSIL)
- High-grade squamous intraepithelial lesions (HSILs)
- Cervical Carcinoma
- Colposcopic or cervicagraphic findings of cervical cancer

Before using the Cervical Rotating Biopsy Punch, clinicians should always consider the need to obtain a tissue diagnosis relative to the risk of contraindications, such as coagulation disorder, severe heart failure, and hypovolemic shock and extreme care should be taken in the presence of active infections such as gonorrhoea, herpes, chlamydia, pelvic inflammatory disease, excessive bleeding other risks include pregnancy, postmenopausal bleeding, anovulatory bleeding, intermenstrual bleeding, or bleeding from a complication of pregnancy.

PATIENT POPULATION

The patient population are females, the target group is to be determined by a suitably qualified clinician or healthcare practitioner.

WARNINGS AND PRECAUTIONS

The device is to be used by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.

- This product is supplied sterile in an undamaged, unopened package as a single use device. Do not resterilise or reuse.
- Before using, inspect package for punctures, tears, security of the seal or any other evidence of the sealed pouch having been compromised prior to placement of the contents in the sterile field.
- DO NOT USE if package is damaged.
- If the device accidentally becomes contaminated before treatment, do not use and dispose accordingly.
- Do not use after the expiration date.
- Make sure the device is compatible with the speculum or any other employed instrument, before use.
- Do not use the device outside the intended use.
- Do not alter the instrument structurally, force it or bend it.
- Beware of the biting jaw being open during introduction into, advancement through and removal from the patients. When the jaws are open, damage to the patient, device and speculum may occur.
- The use of the Cervical Rotating Biopsy Punch on patients can cause varying degrees of pain.
- Instruments made from Stainless Steel may contain Nickel. Although the contact between instruments and patient is transient, it may cause an allergic reaction in patients with sensitivity to this material. Users should risk assess the use of the instruments in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.

POSSIBLE ADVERSE EFFECTS

The use of the Cervical Rotating Biopsy Punch on patients can cause varying degrees of pain. Patient factors might also influence pain. Gynaecologists are advised to prepare and counsel patients prior the procedure that they will experience some form of pain and/or discomfort. Preparing the patient before the procedure may help minimize discomfort.

SHELF LIFE AND STERILITY

Instruments are supplied in a sterile condition, packaged individually and double-wrapped. The instruments are sterilized by exposure to gamma irradiation at a minimum dose of 25 kGy – indicated by the symbol on the label. They remain sterile as long as the package integrity has not been violated.

Packaging must be inspected before use. Do not use any component from an opened or damaged package.

The expiration date is printed on the label. Do not use instruments after the expiration date.

STORAGE AND HANDLING

Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes.

Care must be exercised in handling of wrapped cases or individual instruments to prevent damage to the sterile barrier.

Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

OPERATIONAL USE

Introduce the Cervical Rotating Biopsy Punch through the speculum opening and position it proximal to the cervix area of interest. Rotate the device as necessary using the purple rotational element to ensure optimal position for biopsy and patient interaction. Position the device so the cervical area of interest is between the jaws and pull the trigger to close jaws. Hold the trigger to preserve the biopsy sample just taken. Retrieve the removed segment of tissue for biopsy and collect it accordingly.

DEVICE DISPOSAL

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste. Handle and dispose this product in accordance with accepted medical practice and with applicable local, state and national laws and regulations.

The Cervical Rotating Biopsy Punch does not contain:

- Substances of human origin
- Animal tissues
- Dangerous or hazardous substances
- Radioactive material
- Any substance that may be a Medicinal Product as defined in Article I of Directive 2001/83/EC

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use. If you re-use a device, you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients, including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS.
3. Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting.
 - Other materials: May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

All trademarks herein are the property DTR Medical Ltd, or one of its affiliates unless otherwise indicated.

The words “Authorized Representative” and the name and address of the Authorized Representative if the legal manufacturer of the product and on the product label is from outside of the Member States of the EU.

“Swiss Representative” may also be present if the product is not legally manufactured in Switzerland.

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Summit Medical LLC