# DTR.M.059

## Endoscopic Suction Devices

### Instructions for Use

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## DESCRIPTION

The Endoscopic Suction Devices are a range of suction handles used to suction ear wax, debris, bodily fluids and foreign bodies during ENT procedures.

## MATERIALS OF CONSTRUCTION

Product Component	Details
Metal tube	304 Stainless Steel
Head/Handle	ABS Plastic Hub

## INTENDED USE

The Endoscopic Suction range are suction handles which are used during ENT procedures to remove ear wax, debris, bodily fluids and foreign bodies.

## INTENDED USER

The devices are intended to be used by appropriately trained healthcare professionals.

## INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

Devices are indicated for use in ENT procedures.

These devices are made of Stainless Steel, which may cause an allergic reaction in patients with sensitivity to Nickel. User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to Nickel.

## PATIENT POPULATION

Anyone of any age in the human population as determined by a suitably qualified healthcare professional.

## WARNINGS AND PRECAUTIONS

- This device is supplied STERILE and ready to use.

- The device is for **SINGLE USE ONLY**. Do **NOT** re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- If the device accidentally becomes contaminated before use, do not use and dispose of accordingly.
- Care should be taken when using the device to prevent damage to the tissues and structures.

## **SHELF LIFE AND STERILITY**

- The instruments are sterilised by exposure to Ethylene Oxide, as indicated by the symbol on the label. They remain sterile as long as the package integrity has not been compromised.
- The expiration date is printed on the label. Do not use instruments after the expiration date. The shelf life is 5 years.
- Devices are single use, transient devices. Therefore, the device lifetime in use is less than 60 minutes continuous use.

## **STORAGE AND HANDLING**

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes.
- Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

## **TRANSPORTATION**

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

## **OPERATIONAL USE**

This is a standard device, for operational use please refer to indications and intended use.

## **DEVICE DISPOSAL**

After use the device should be disposed of in accordance with the healthcare settings disposal procedure. The product comes into

contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

## **REPORTING PROBLEMS**

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at [www.dtrmedical.com/contact](http://www.dtrmedical.com/contact) and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: [www.tga.gov.au](http://www.tga.gov.au)

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC