DTR.M.063 FESS SINUS SUCTION TUBES

Instructions for Use

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74-2060, 74-2063

DESCRIPTION

The FESS Sinus Suction Tube is used to remove bodily fluids and foreign objects in the nasal cavity.

MATERIALS OF CONSTRUCTION

Product Component	Details
Metal tube	304 Stainless Steel
Head/Handle	ABS Plastic Hub

INTENDED USE

The Functional Endoscopic Sinus Surgery (F.E.S.S) Suction Tube are suction handles which are used during ENT procedures to remove bodily fluids and foreign bodies.

INTENDED USER

The device is intended to be used by appropriately trained surgeons.

INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

FESS Sinus Suction Tube is indicated for use in the nasal cavity for nasal procedures.

These devices are made of Stainless Steel, which may cause an allergic reaction in patients with sensitivity to Nickel. User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to Nickel.

PATIENT POPULATION

Anyone of any age in the human population as determined by a suitably qualified surgeon.

WARNINGS AND PRECAUTIONS

- This device is supplied STERILE and ready to use.

- The device is for **SINGLE USE ONLY**. Do **NOT** re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- If the device accidentally becomes contaminated before use, do not use and dispose of accordingly.
- Care should be taken when using the device to prevent damage to the tissues and structures of the nasal cavity, ear canal and drum.

SHELF LIFE AND STERILITY

- The instruments are sterilised by exposure to Ethylene Oxide, as indicated by the symbol on the label. They remain sterile as long as the package integrity has not been violated.
- Packaging must be inspected before use. Do not use any component from an opened or damaged package.
- The expiration date is printed on the label. Do not use instruments after the expiration date. The shelf life is 5 years.
- Devices are single use, transient devices. Therefore, the device lifetime in use is less than 60 minutes continuous use.

STORAGE AND HANDLING

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes. Store in a cool, dry place.
- Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

POSSIBLE ADVERSE EFFECTS

No known adverse events from use of the device as indicated

OPERATIONAL USE

This is a standard device, for operational use please refer to indications and intended use.

DEVICE DISPOSAL

After use the device should be disposed of in accordance with the healthcare settings disposal procedure. The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC