DTR.M.067

NON STERILE SILICONE TUBING

Instructions for Use

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DESCRIPTION

The Non Sterile Silicone Tubing is intended to be cut to the desired length for use outside of the sterile surgical field.

MATERIALS OF CONSTRUCTION

Product Component	Details
Tubing	Silicone

INTENDED USE

The Non Sterile Silicone Tubing is intended to be cut to the desired length for use outside of the sterile surgical field.

INTENDED USER

The devices are intended to be used by appropriately trained healthcare professionals.

INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

The Non Sterile Silicone Tubing is indicated to be cut to the desired length for use outside of the sterile surgical field.

These devices are made of silicone which may cause an allergic reaction in users with sensitivity to this material. User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with users with known sensitivity to silicone.

PATIENT POPULATION

Not applicable as these devices are non-patient contacting.

WARNINGS AND PRECAUTIONS

- The device is supplied NON-STERILE. Do not use in the sterile surgical field.
- The device is for SINGLE USE ONLY, do NOT re-use.
- The device is non patient contacting.
- Do not use if the packaging has been opened or damaged.

- If the device accidentally becomes contaminated before use, do not use and dispose of accordingly.

SHELF LIFE AND STERILITY

These devices are not intended to be used on a patient. The devices are non-patient contacting and continuous time in use is less than 24 hours. The expiration date is printed on the label. Do not use instruments after the expiration date. The shelf life is 5 years.

STORAGE AND HANDLING

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

OPERATIONAL USE

This is a standard device, for operational use please refer to indications and intended use.

DEVICE DISPOSAL

After use the device should be disposed of in accordance with the healthcare settings disposal procedure. The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has

occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

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DTR Medical Ltd, 17 Clarion Court,