# DTR.M.072

## Micro Suction Devices/Aspiration Cannulae

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## DESCRIPTION

The Micro Suction Handles / Aspiration Cannula are used within general, orthopaedic, maxillofacial and ENT procedures to remove bodily fluids, ear wax and foreign bodies.

## MATERIALS OF CONSTRUCTION

Product Component	Details
Metal tube	304 Stainless Steel
Head/Handle: 74-1260; 74-1262A; 74-1230; 74-1231; 74-1231A; 74-1232; 74-1232A; 74-1233; 74-1234; 74-1236; 74-1264; 74-1265;	Polypropylene Hub
Head/Handle: 74-1268; 74-1297	PVC
Head/Handle: 74-1267;	PC

## INTENDED USE

The Micro Suction Devices / Aspiration Cannulae are intended for use in General, Orthopaedic, Maxillofacial and ENT procedures.

## INTENDED USER

The device is intended to be used by appropriately trained healthcare professionals.

## INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

The Micro Suction Devices / Aspiration Cannulae are indicated for the removal of bodily fluids, ear wax and foreign bodies.

Do not use with in direct contact with the heart, central circulatory (arteriae pulmonales, aorta ascendens, arcus aortae, aorta descendens to the bifurcatio aortae, arteriae coronariae, ateria carotis communis,

arteria carotis externa, arteria carotis interna, arteriae cerebrales, truncus brachiocephalicus, venae cordis, venae pulmonales, vena cava superior and vena cava inferior) and central nervous system (brain, meninges and spinal cord).

These devices are made of Stainless Steel, which may cause an allergic reaction in patients with sensitivity to Nickel. User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to Nickel.

These devices are not indicated to be used for irrigation.

## **PATIENT POPULATION**

Anyone of any age in the human population as determined by a suitably qualified healthcare professional.

## **WARNINGS AND PRECAUTIONS**

- This device is supplied STERILE and ready to use.
- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- If the device accidentally becomes contaminated before use, do not use and dispose of accordingly.
- Care should be taken when using the device to prevent damage to the tissues and structures.

## **SHELF LIFE AND STERILITY**

- The instruments are sterilised by exposure to Ethylene Oxide, as indicated by the symbol on the label. They remain sterile as long as the package integrity has not been compromised.
- The expiration date is printed on the label. Do not use instruments after the expiration date. The shelf life is 5 years.
- Devices are single use, transient devices. Therefore, the device lifetime in use is less than 60 minutes continuous use.

## **STORAGE AND HANDLING**

- Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes.

- Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.
- Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

## **TRANSPORTATION**

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

## **OPERATIONAL USE**

This is a standard device, for operational use please refer to indications and intended use.

## **DEVICE DISPOSAL**

After use the device should be disposed of in accordance with the healthcare settings disposal procedure. The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

## **REPORTING PROBLEMS**

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at [www.dtrmedical.com/contact](http://www.dtrmedical.com/contact) and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred. For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: [www.tga.gov.au](http://www.tga.gov.au).

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC