**DTR.M.015**  
**CROCODILE**  
**MICRO FORCEPS,**  
**CROCODILE FINE JAW**  
**MICRO FORCEPS &**  
**CROCODILE FORCEPS**  
Instructions for Use

---

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**DESCRIPTION**

Crocodile Forceps are sterile single-use stainless steel devices used for grasping and holding. They have a rigid body and ring handles. The body is thin with a 117° bend where it tapers to a tip. At the tip is a pair of serrated jaws which are used to grasp and hold. DTR Medical sells a range of sizes to fit the ENT procedures. The device is sold in a normal steel or a matte black finish. The matte black finish prevents reflection from illumination devices used in ENT procedures.

**MATERIALS OF CONSTRUCTION**

Product Component	Details
Crocodile Forceps	410 Stainless Steel

**INTENDED USE**

The forceps are used in multiple surgical ENT procedures for both grasping and holding tissue in difficult to reach areas.

**INTENDED USER**

A health care professional with the required training to use the forceps.

## **INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS**

Before using the forceps ensure the Tympanic membrane is intact.

## **PATIENT POPULATION**

Anyone of any gender or any age in the human population that is identified as needing an ENT procedure.

## **WARNINGS AND PRECAUTIONS**

The device is to be used by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.

- This product is supplied sterile in an undamaged, unopened package as a single use device. Do not resterilise or reuse.
- Before using, inspect package for punctures, tears, security of the seal or any other evidence of the sealed pouch having been compromised prior to placement of the contents in the sterile field.
- **DO NOT USE** if package is damaged.
- If the device accidentally becomes contaminated before treatment, do not use and dispose accordingly.
- Do not use after the expiration date.
- Make sure the device is compatible with the speculum or any other employed instrument, before use.
- Do not use the device outside the intended use.
- Do not alter the instrument structurally, force it or bend it.
- Beware that inappropriate use can cause tissue damage.
- Instruments made from Stainless Steel may contain Nickel. Although the contact between instruments and patient is transient, it may cause an allergic reaction in patients with sensitivity to this material. Users should risk assess the use of the instruments in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- Ensure all packaging is removed before use, including tip shields or caps on the product.

## **SHELF LIFE AND STERILITY**

Instruments are supplied in a sterile condition, packaged individually and single-wrapped. The instruments are sterilized by exposure to

Ethylene Oxide (EtO) sterilisation process.

Packaging must be inspected before use. Do not use any component from an opened or damaged package. The expiration date is printed on the label. Do not use instruments after the expiration date.

## **STORAGE AND HANDLING**

Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes. Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.

Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

## **TRANSPORTATION**

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

## **OPERATIONAL USE**

Crocodile forceps use finger loops to facilitate handling. The finger loops are usually grasped by the thumb and middle or ring fingers, while the index finger helps stabilise and guide the instrument. Once the instrument is in place, pressure can be applied to close the jaws.

## **DEVICE DISPOSAL**

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste. Handle and dispose this product in accordance with accepted medical practice and with applicable local, state and national laws and regulations.

The Crocodile Micro Forceps, Crocodile Fine Jaw Micro Forceps and Crocodile Forceps does not contain:

- Substances of human origin
- Animal tissues
- Dangerous or hazardous substances

- Radioactive material
- Any substance that may be a Medicinal Product as defined in Article I of Directive 2001/83/EC

## **REPORTING PROBLEMS**

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at [www.dtrmedical.com/contact](http://www.dtrmedical.com/contact) and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: [www.tga.gov.au](http://www.tga.gov.au).

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC