DTR.M.017 ORMEROD FORCEPS

Instructions for Use

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DESCRIPTION

Ormerod Forceps are sterile single-use stainless steel devices. They have a rigid body with ring handles. The body is thin with a 117° bend where it tapers to a tip. At the tip is a pair of cup shaped jaws with a small cavity within, which are used to grasp and hold.

MATERIALS OF CONSTRUCTION

Product Component	Details
Ormerod Forceps	410 Stainless Steel

INTENDED USE

The forceps are used in multiple surgical ENT procedures for both grasping and holding tissue in difficult to reach areas.

INTENDED USER

A health care professional with the required training to use the forceps.

INDICATIONS, CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

Before using the forceps ensure the Tympanic membrane is intact.

PATIENT POPULATION

Anyone of any gender or any age in the human population that is identified as needing an ENT procedure.

WARNINGS AND PRECAUTIONS

The device is to be used by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.

- This product is supplied sterile in an undamaged, unopened package as a single use device. Do not resterilise or reuse.
- Before using, inspect package for punctures, tears, security of the seal or any other evidence of the sealed pouch having been compromised

- prior to placement of the contents in the sterile field.
- DO NOT USE if package is damaged.
 - If the device accidentally becomes contaminated before treatment, do not use and dispose accordingly.
 - Do not use after the expiration date.
 - Make sure the device is compatible with the speculum or any other employed instrument, before use.
 - Do not use the device outside the intended use.
 - Do not alter the instrument structurally, force it or bend it.
 - Beware that inappropriate use can cause tissue damage.
 - Instruments made from Stainless Steel may contain Nickel. Although the contact between instruments and patient is transient, it may cause an allergic reaction in patients with sensitivity to this material. Users should risk assess the use of the instruments in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
 - Ensure all packaging is removed before use, including tip shields on the product.

SHELF LIFE AND STERILITY

Instruments are supplied in a sterile condition, packaged individually and single-wrapped. The instruments are sterilized by exposure to Ethylene Oxide (EtO) sterilisation process.

Packaging must be inspected before use. Do not use any component from an opened or damaged package. The expiration date is printed on the label. Do not use instruments after the expiration date.

STORAGE AND HANDLING

Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes. Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.

Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This

packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

OPERATIONAL USE

Ormerod forceps use finger loops to facilitate handling. The finger loops are usually grasped by the thumb and middle or ring fingers, while the index finger helps stabilise and guide the instrument. Once the instrument is in place, pressure can be applied to close the jaws.

DEVICE DISPOSAL

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste. Handle and dispose this product in accordance with accepted medical practice and with applicable local, state and national laws and regulations.

The Ormerod Forceps do not contain:

- Substances of human origin
- Animal tissues
- Dangerous or hazardous substances
- Radioactive material
- Any substance that may be a Medicinal Product as defined in Article I of Directive 2001/83/EC

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC