DTR.M.025 ARTERY FORCEPS

Instructions for Use

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
Baby Overholt
Artery Forceps

REF

CFS2001-D	Crile Artery Forceps, 145mm Straight
CFC2001-D	Crile Artery Forceps, 145mm Curved on Flat
KAF2014-D	Kocher Artery Forceps 140mm
KOF2181	Kocher-Ochsner Forceps, Straight 185mm
MQC2012	Micro Mosquito Forceps 120mm Curved
MQT2001	Mosquito Forceps
MQC2025-D	Mosquito Forceps Curved 125mm
MQC2125-D	Mosquito Forceps Curved Toothed 1x2 125mm
SWF2010	Spencer Wells Artery Forceps 130mm straight
SWF2015	Spencer Wells Artery Forceps 155mm straight
SWF2018	Spencer Wells Artery Forceps 185mm straight
SWF2020	Spencer Wells Artery Forceps 200mm straight
BOC2014	Baby Overholt Artery Forceps Curved, 140mm

**DESCRIPTION**

The Artery forceps are hinged devices that are designed for grasping vessels and temporarily suspending blood flow to allow a surgical procedure to be conducted.

MATERIALS OF CONSTRUCTION

Product Component	Details
BOC2014	420 Stainless Steel
CFC2001-D	410 Stainless Steel
CFS2001-D	410 Stainless Steel
KAF2014-D	410 Stainless Steel
KOF2181	410 Stainless Steel
MQC2012	420 Stainless Steel
MQC2025-D	410 Stainless Steel

MQC2125-D	420 Stainless Steel
MQT2001	420 Stainless Steel
SWF2010	410 Stainless Steel
SWF2015	410 Stainless Steel
SWF2018	410 Stainless Steel
SWF2020	410 Stainless Steel

INTENDED USE/INDICATIONS FOR USE

The Artery forceps are used to grasp vessels and apply sufficient pressure to temporarily suspend blood flow to allow a surgical procedure to be conducted.

INTENDED USER

The device is intended to be used by appropriately trained healthcare professionals.

CONTRAINDICATIONS AND PATIENT SELECTION FACTORS

These devices are made of Stainless Steel, which may cause an allergic reaction in patients with sensitivity to Nickel. The User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to Nickel.

PATIENT POPULATION

Anyone in the human population as determined by a suitable qualified healthcare professional.

WARNINGS AND PRECAUTIONS

The device is to be used by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.

- This product is supplied sterile in an undamaged, unopened package as a single use device. Do not resterilise or reuse.
- Before using, inspect package for punctures, tears, security of the seal or any other evidence of the sealed pouch having been compromised

- prior to placement of the contents in the sterile field.
- **DO NOT USE** if package is damaged.
- If the device accidentally becomes contaminated before treatment, do not use and dispose accordingly.
- Do not use after the expiration date.
- Do not use the device outside the intended use.
- Do not alter the instrument structurally, force it or bend it.
- Beware that inappropriate use can cause tissue damage.

SHELF LIFE AND STERILITY

Instruments are supplied in a sterile condition, packaged individually. The instruments are sterilized by exposure to Ethylene Oxide (EtO) sterilisation process.

Packaging must be inspected before use. Do not use any component from an opened or damaged package. The expiration date is printed on the label. Do not use instruments after the expiration date.

STORAGE AND HANDLING

Sterile, packaged instruments should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, sunlight and temperature/humidity extremes. Care must be exercised in handling of individual instruments to prevent damage to the sterile barrier.

Sterile instrument packages should be carefully examined prior to opening to ensure that package integrity has not been compromised.

TRANSPORTATION

The instruments are supplied in their designated packaging. This packaging ensures that every instrument is kept in a way that they do not receive any damage and that their functionality is preserved during transportation.

OPERATIONAL USE

Artery forceps use finger loops to facilitate handling. The finger loops are usually grasped by the thumb and middle or ring fingers, while the index finger helps stabilise and guide the instrument. Once the instrument is in place, pressure can be applied to close the jaws.

DEVICE DISPOSAL

Used devices should be disposed of as clinical waste. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste. Handle and dispose this product in accordance with accepted medical practice and with applicable local, state and national laws and regulations.

The Artery Forceps do not contain:

- Substances of human origin
- Animal tissues
- Dangerous or hazardous substances
- Radioactive material
- Any substance that may be a Medicinal Product as defined in Article I of Directive 2001/83/EC

REPORTING PROBLEMS

Healthcare professionals, users and patients should report any suspected serious incident related to the device by informing the manufacturer at www.dtrmedical.com/contact and the competent authority, ministry of health, or delegated agency in which the suspected serious incident has occurred.

For patients in Australia please visit the Therapeutic Goods Administration (TGA) website: www.tga.gov.au.

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DTR Medical Ltd, 17 Clarion Court,
Summit Medical LLC