



NEW ENT Products 2020




Lempert Suction Handle

140mm

- 142° angled handle for improved visibility
- Lightweight device for ease of use
- Rotation for improved visibility
- Thumb control for suction
- Single packed for fast access

LSH2001-S



Thomassin Dissector

180mm

- Used for removing tissue during Endoscopic Ear Surgery
- Double ended for increased flexibility
- Angled to allow tips to operate on deep and difficult to reach areas in the middle ear
- Thin shaft provides clear view of the surgical site
- Diamond knurled handle provides grip during the procedure

TDD1018




Wire Cutting Scissors

125mm

- New and improved Scissors enable the cutting of a wide variety of wire
- Used in jaw fixation
- Allows patients or relatives to cut the wire if any vomiting or choking are a risk
- Ideal for cutting standard Cerclage wire efficiently along the whole length of the blade

WCSH2001





Endoscopic House Curette Steep Angle




Endoscopic House Curette

180mm

- Engineered with hardened steel stronger than titanium
- Provides a better cutting edge with increased durability
- Ensures less pressure is required when in use
- 3 sharp, angled scoop sizes to allow for a clean cut
- Clearly labelled 2.0 & 3.0mm ends
- Used for curettage of bone and soft tissue within the middle ear
- Slender shaft which does not obstruct the user's view of the surgical site

HCDL2230, HCDM2230, HCDS2230



Keele Mucosectomy Forceps

250mm

- Ideal for use during Trans Oral Mucosectomy procedure
- Forceps are used to hold, control, manipulate and stabilise tissue in the oral cavity and manoeuvre the tongue during cauterisation
- 140° angled tip for enhanced access
- 1x2 Teeth to maintain firm specimen retention

KMF1025

SuperCut Iris Scissors

90mm & 110mm

- One serrated and one sharp cutting blade allows a better grip and cut than traditional alternatives
- Sharp tips make the Scissors ideal for precise cutting of soft tissue, sutures and wound dressings
- Available with straight or curved blades
- Black handles make the instrument distinct

SISS2111, SISC2111, SISS2110, SISC2110



SuperCut Iris Scissors 90mm



YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

