


Innovia Face Visor Specification Sheet



Product Overview

- Designed to provide face protection from aerosols and minimise airborne cross-contamination
- CE marked to regulation 2016/425 (Notified Body BSI Netherlands Group BV)
- Meets all applicable requirements of EN 166:2002 Personal eye Protection
- Provides top, side, and front face protection for the healthcare professional
- 22.8cm long screen provides maximum protection.
- Elastic headband provides secure retention and easy application
- Can be worn over prescription eyewear and facemasks
- Anti-fog and optically clear screen for maximum visibility
- Lightweight and comfortable
- Latex-free
- Manufactured in an ISO Class 9 or higher classification UK cleanroom
- Non-sterile, single use, individually wrapped

Description	Pack Size	Box Size	Code
Innovia Face Visor	1x Visor	100x Packs	CM-5001

Innovia Face Visor Specification Sheet

Material Specification

Product Component

Specification

Supportive Headrest	Adhesive Backed RH-32 Polyurethane Foam
Visor	PET Melinex, Polyester, Clear
Headband	Latex Free Elastic
Individual Face Visor Packaging/ Bag	Grip Seal, Medium Square (12.75" x 12.75")
Outer Box/Carton Packaging	Plain Brown, Single Wall (610x457x457)

Intended Use

The device is intended to provide top, side, and front face protection from aerosols and to minimize airborne cross-contamination.

This protective Face Visor is manufactured for COVID-protection only. This protective Face Visor is not a PPE device for general use and shall not be used for purposes other than protection against COVID-19.

Instructions for Use

Instructions for use are provided.

Conformity to the European Regulations

This device is Category II PPE in accordance with regulation (EU) 2016/425 of the European Parliament and of the Council.