



Semi-sharp=0.7mm tip



Secure Hold Without the Trauma with our Gillies Skin Hook



Sharp=narrows to a
fine and Sharp point



The Design Checklist

- ✓ Available in both Sharp and Semi-sharp the 182mm Gillies Skin Hook offers you flexibility in a variety of procedures
- ✓ The Sharp Hook enables enhanced accuracy, grip and precision
- ✓ The Semi-sharp Hook is suitable for holding skin whilst taking a biopsy, minimising patient trauma
- ✓ Manufactured using high-quality stainless steel to ensure first-time sharpness
- ✓ Our Skin Hooks are 100% quality inspected to give you confidence

Clinicians' Experiences

Don't just take our word for it ...

"Amazing blunt Hook which I have found perfect for safely isolating and gently moving small blood vessels as well as using as a gentle Semi-sharp Gillies Skin Hook. Very well made and lovely to handle."

Prof. Mark Whiteley MS FRCS(Gen) FCPHleb, Consultant Venous Surgeon and Consultant Phlebologist, The Whiteley Clinic

Clinical Value - Sharp & Semi-sharp to Improve Patient Experience



TIME

With prompt delivery guaranteed, we endeavour to meet your service needs. Avoiding costly delays, reducing hospital stock holding and improving lead time for instruments



LIFE

The Sharp Hook removes the risk of fine tips being blunted through repeated sterilisation, meaning a precise retraction

The Semi-sharp Hook ensures increased comfort as the atraumatic tips mean reduced tissue damage and less patient pain



COST

Our Single-use Gillies Skin Hooks are reassuringly Sharp to enable you to achieve first time results. This is a cost-effective alternative in comparison to poor quality, cheap instruments that require opening multiple packs per patient

No burden of Sterilisation costs