

## Reusable Innovia Face Visor Product Specification



### The Design Checklist

- Manufactured to highest PPE regulation requirements**  
 CE Category III 2016/425 as per Annex II (Notified Body Satra Technology Europe Ltd)
- Reusable design reduces environmental impact**  
 Allows for up to 5 days of safe wear by the same individual
- N rated hydrophilic anti-fog coating**  
 Prevents condensation forming for an optically clear screen and maximum visibility (compliant with EN166:2001 Clause 7.3.2)
- Permanent marker compatible Visor includes space for individual name and customisation**
- Pre-printed 5-day tally tick box feature allows wearer to keep track of number of days used and when to dispose of and replace**
- Ultra-soft grooved foam headrest for increased comfort and removal of dermal heat**
- Designed to provide face protection from aerosols and minimise airborne cross-contamination**
- Meets all applicable requirements of EN 166:2001 Personal Eye Protection**
- Provides top, side, and front face protection for the healthcare professional**
- 22.8cm long screen provides maximum protection**
- Elastic headband provides secure retention and easy application**
- Can be worn over prescription eyewear and facemasks**
- Anti-glare performance for optimum clarity within clinical environment**
- Lightweight and comfortable**
- Latex-free**
- Manufactured in an ISO Class 9 or higher UK cleanroom**
- Non-sterile, individually wrapped**

## Reusable Innovia Face Visor Product Specification

### Intended Use

The device is intended to provide top, side, and front face protection from aerosols and to minimise airborne cross-contamination.

This protective Face Visor is manufactured for COVID-protection only. This protective Face Visor is intended to provide protection against liquid splashes.

This Visor is not resistant against spark, flame, or impact.

### Instructions For Use (IFU)

Instructions for use are provided.

### Conformity to European Regulations

This device is Category III PPE in accordance with regulation (EU) 2016/425 of the European Parliament and of the Council.

This device complies with the requirements of EN 166:2001 and associated standards and bear the following marks:

Marking	Description
1	Optical Class
3	Field of use: Liquid Product protects against liquid droplets or splashes
N	Resistance to fogging



### Product Specification

Description	Pack Size	Box Size	Code
Reusable Innovia Face Visor	1x Visor	100x	CM-5001R

### Material Breakdown

Product Component	Specification
Supportive Headrest	Adhesive Backed RH-32 Polyurethane Foam
Visor	250µm Polycarbonate, Clear
Headband	Latex Free Elastic
Individual Face Visor Packaging/ Bag	Grip Seal, Medium Square (12.75" x 12.75")
Outer Box/Carton Packaging	Plain Brown, Double Wall (610x457x457)