




DTR Medical®



AN
INNOVIA
MEDICAL
COMPANY

New ENT Products 2021

Endoscopic House Curette



Endoscopic House Curette Steep Angle



- Engineered with hardened steel stronger than titanium
- Provides a better cutting edge with increased durability
- Ensures less pressure is required when in use
- 3 sharp, angled scoop sizes to allow for a clean cut
- Clearly labelled 2.0 & 3.0mm ends
- Used for curettage of bone and soft tissue within the middle ear
- Slender shaft which does not obstruct the user's view of the surgical site

180mm

Keele Mucosectomy Forceps



- Ideal for use during Trans Oral Mucosectomy procedure
- Forceps are used to hold, control, manipulate and stabilise tissue in the Oral Cavity and manoeuvre the tongue during Cauterisation
- 140° angled tip for enhanced access
- 1x2 Teeth to maintain firm specimen retention

250mm

TO EXPERIENCE THE PRECISION OF DTR MEDICAL STERILE SINGLE-USE INSTRUMENTS PLEASE CONTACT US

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com



SuperCut Iris Scissors

- One serrated and one sharp cutting blade allows a better grip and cut than traditional alternatives
- Sharp tips make the Scissors ideal for precise cutting of soft tissue, sutures, and wound dressings
- Available with straight or curved blades
- Black handles make the instrument distinct

90mm & 110mm



SuperCut Iris Scissors 90mm



SNAP, Endoscope Guide

- Protects from airborne risks and diseases
- Designed with the aim to keep Clinicians safe and so better able to deliver high quality care to their patients
- Facilitates the movement of a Nasendoscope up to 4mm in diameter
- Ensures safe passage of a Nasendoscopy through a Surgical Mask, creating an aperture for the Endoscope to pass through into the Nasal Cavity
- Scientifically proven to reduce the prevalence of AGP (Aerosol Generating Particles) during examination
- Provides better safety and care for both the Patient & Clinician



TO EXPERIENCE THE PRECISION OF DTR MEDICAL STERILE SINGLE-USE INSTRUMENTS PLEASE CONTACT US

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

