



DTR Medical®

AN
INNOVIA
MEDICAL
COMPANY


Keele Mucosectomy Forceps

Carcinomas at the back of the tongue can be removed using one of three procedures:

- TORS – Transoral Robotic Surgery
- TOLS – Transoral Laser Surgery
- TOES – Transoral Endoscopic Surgery

Limitations:

- TORS – Limited view and expensive
- TOLS – Narrow view when using a laryngoscope
- TOES – Carcinomas at the back of the tongue are not easily accessible due to lack of suitably sized instruments



Keele Mucosectomy Forceps are used in TOES:

- The 250mm length and 140° angled tip provides the user with extended reach to grasp carcinomas behind the tongue
- 1x2 teeth to manoeuvre and manipulate tissue
- Used alongside an electrode with a ceramic malleable tip which can be bent to the same angle as the forceps
- Quick set up in comparison to robotic surgery which can double the length of an operation
- Reduce cost to perform procedure
- Wider view of the oral cavity

YOUR SINGLE-USE SOLUTION WITH CONFIDENCE

t +44 (0)1792 797910 e marketing@dtrmedical.com w www.dtrmedical.com

