


## 'Endo' Your Search Here for a High-Quality Endometrial Biopsy with our Sterile Single-use Medgyn Endosampler



### The Design Checklist

- ✓ Designed to facilitate sampling for Gynaecological procedures such as Cytologic and Histologic diagnosis
- ✓ Enables Endometrial scraping and aspiration
- ✓ Conveniently bent with a sharp slot on the side
- ✓ Semi-ridged 3mm curved Curette and 10CC Syringe
- ✓ Features luer lock spring mechanism for larger and multiple sampling
- ✓ Engineered to prevent any potential backflow of the collected biopsy



See our video demonstration of the Medgyn Endosampler at:  
[www.dtrmedical.com/products/medgyn-products/](http://www.dtrmedical.com/products/medgyn-products/)

Or scan  
this code



## Clinical Value - Struggling for an Accurate First Time Sample?



### TIME

Reliable supply, supported by strong service



### LIFE

Allows you to take an accurate sample from all areas of the Endometrium

Reliability when taking an Endometrial tissue sample

A safe procedure, well tolerated by patients with a low risk of perforation

Avoids sample contamination



### COST

Achieve superior first-time accuracy

No burden of sterilisation costs

## Clinicians' Experiences - Don't just take our word for it ...

**"The Endosampler is not only an effective kit for outpatient Endometrial Biopsy, it also has added benefits:**

- 1. Occasionally small Polyps can be dislodged during Endometrial Biopsy**
- 2. Occasionally help during retrieval of retained Coil when threads are in the Uterus**
- 3. Can be used to evacuate saline from Uterine Cavity and Endometrial sampling just before Novasure Endometrial ablation, 1 instrument for 2 purposes in this situation."**

Gynaecology Department, Calderdale and Huddersfield NHS Foundation Trust