

**Extended Range to Suit your Colposcopy Requirements -
Try our Single-use DTR Medical[®] LLETZ Colposcopy Loops & Balls**



Colposcopy LLETZ Loop Electrodes

Colposcopy LLETZ Ball Electrodes

The Design Checklist

Colposcopy LLETZ Loop Electrodes

- ✓ Tungsten Loops provide ideal tension for greater control
- ✓ Stiffer at smaller dimensions which prevents the wire from folding when cutting
- ✓ Universal shaft fits most electrosurgical pencils

Colposcopy LLETZ Ball Electrodes

- ✓ Gap between handle and Ball Electrode for better access
- ✓ Two sizes available in 130mm & 152mm length, both with the option of 3mm or 5mm Ball diameter
- ✓ 152mm variants suitable for larger patient
- ✓ Non-stick Teflon coating changes the heat dissipation and prevents tissue from sticking to the instrument

Clinical Value – Assisting your Busy Colposcopy Clinic



TIME

Recognised colour coding with Loop size printed on the handle to easily distinguish between the variants

With prompt delivery guaranteed, we endeavour to meet your service needs. Avoiding costly delays, reducing hospital stock holding and improving lead time for Instruments

No time wasted waiting for sterilisation to return, avoiding possible clinic delays



LIFE

High cutting performance with reduced thermal damage to tissue

Single-use ensures first time sharpness and performance with no cross infection risks or health and safety standards



COST

No burden of sterilisation costs



LLETZ Loop Electrodes Handle Size Recognition Guide			
Size (mm)	Length (mm)	Handle	Loop Colour Code
5 x 5	130	Purple	Purple
10 x 5	130	Light Purple	Light Purple
10 x 10	130	Yellow	Yellow
15 x 8	130	Grey	Grey
15 x 12	130	Green	Green
15 x 15	130	Orange	Orange
20 x 10	130	Red	Red
20 x 12	130	White	White
20 x 15	130	Light Blue	Light Blue
20 x 20	130	Pink	Pink
25 x 8	130	Tan	Tan