

The 'Pipette' of the Bunch – Use our Sterile Single-use Medgyn Pipette IV Hole



The Design Checklist

- ✓ Suitable for a range of Gynaecological procedures including Cytologic and Histologic diagnosis
- ✓ Designed with enhanced flexibility and longitudinal strength
- ✓ Semi-rigid, tapered 3.1mm cannula enables a smooth insertion and minimises the possibility of perforation
- ✓ Luer lock when multiple biopsies are required

Clinical Value – Superior Sampling with our Pipette IV Hole



TIME

Reliable supply, supported by strong service



LIFE

Improved sampling when used to clearly differentiate Endometrial Tissue without the need for Anaesthesia in an outpatient setting

Flexibility assures better sampling from all areas of the Uterine Cavity, even from difficult-to-reach cornual areas

Smooth end reduces patient discomfort

Avoids sample contamination



COST

Achieve first-time accuracy

No burden of sterilisation costs

Clinicians' Experiences - Don't just take our word for it, contact us to trial the Medgyn Pipette IV Hole in your Hospital