

Denver Splint[®]

Nasal Splint




KEY PRODUCT FEATURES

- ✓ Stabilization of the nose after surgery
- ✓ Reduction of edema after nasal surgery or trauma
- ✓ Moldable aluminum retains nasal bones in position
- ✓ Layers give cushioned pressure on subcutaneous tissue
- ✓ Dorsal pad reduces dead space and possibility of polybeak
- ✓ Approximately 60 seconds to apply
- ✓ Adhesive
- ✓ No water mess
- ✓ Inconspicuous
- ✓ Water resistant
- ✓ 4 Sizes, One size per box
- ✓ Kit and Non-Kit

KIT INCLUDES

- One Denver Splint[®]
- Two Pre-cut Protecto Strips, sterile
- One alcohol wipe, sterile
- One skin prep, sterile
- One dorsal pad



NON-KIT INCLUDES

- One Denver Splint[®]
- One dorsal pad

Denver Splint®

Nasal Splint

Denver Splint® - Two Piece

Denver Splint® Series 1500

Soft moldable aluminum and Velcro™

Code	Description	Dimensions mm	QTY/Box
10-1500-05KM	Mini	50 x 28	Kit 5
10-1500-05KP	Petite	58 x 30	Kit 5
10-1500-05KS	Small/ Medium	65 x 35	Kit 5
10-1500-05KL	Large	69 x 43	Kit 5
10-1500-10NKM	Mini	50 x 28	Non-kit 10
10-1500-10NKP	Petite	58 x 30	Non-kit 10
10-1500-10KNS	Small/ Medium	65 x 35	Non-kit 10
10-1500-10NKL	Large	69 x 43	Non-kit 10

Denver Splint® Series 5000

Soft moldable aluminum and Velcro™ with attached dorsal pad

Code	Description	Dimensions mm	QTY/Box
10-5000-05KM	Mini	50 x 28	Kit 5
10-5000-05KP	Petite	58 x 30	Kit 5
10-5000-05KS	Small/ Medium	65 x 35	Kit 5
10-5000-05KL	Large	69 x 43	Kit 5
10-5000-10NKM	Mini	50 x 28	Non-kit 10
10-5000-10NKP	Petite	58 x 30	Non-kit 10
10-5000-10KNS	Small/ Medium	65 x 35	Non-kit 10
10-5000-10NKL	Large	69 x 43	Non-kit 10


Denver Splint® - One Piece

Soft moldable aluminum and 1/8" ventilating foam

Denver Splint® Series 2000

Soft moldable aluminum and Velcro™

Code	Description	Dimensions mm	QTY/Box
10-2000-05KM	Mini	50 x 28	Kit 5
10-2000-05KP	Petite	58 x 30	Kit 5
10-2000-05KS	Small/ Medium	65 x 35	Kit 5
10-2000-05KL	Large	69 x 43	Kit 5
10-2000-10NKM	Mini	50 x 28	Non-kit 10
10-2000-10NKP	Petite	58 x 30	Non-kit 10
10-2000-10KNS	Small/ Medium	65 x 35	Non-kit 10
10-2000-10NKL	Large	69 x 43	Non-kit 10



Denver Splint®

Thermoplastic Nasal Splint




KEY PRODUCT FEATURES


- ✓ Rigid, opaque, lightweight, low temperature thermoplastic splinting material
- ✓ Simply soften with hot water and mold over taped nose
- ✓ Stabilization of the nose after surgery
- ✓ Reduction of edema after nasal surgery or trauma
- ✓ Length of wear 7-10 days
- ✓ 4 Sizes, One size per box
- ✓ Available in:
 - Pre-cut in kit and Non-Kit
 - 3×3" Squares
 - 9×12" Sheets (Cut custom shape)

KIT INCLUDES

- One Denver Splint®
- Two Pre-cut Protecto Strips, sterile
- One alcohol wipe, sterile
- One skin prep, sterile
- One dorsal pad



INDICATE SIZE



Denver Splint®

Thermoplastic Nasal Splint

Thermoplastic Denver Splint®				
Pre-cut Denver Splint® Series 4000				
Code	Description	Dimensions mm	Color	QTY/Box
10-4021-05KM	Mini	50 x 28	Ivory 1/16" (1.6mm)	Kit 5
10-4021-05KP	Petite	58 x 30	Ivory 1/16" (1.6mm)	Kit 5
10-4021-05KS	Small	65 x 35	Ivory 1/16" (1.6mm)	Kit 5
10-4021-05KL	Large	69 x 43	Ivory 1/16" (1.6mm)	Kit 5
10-4021-10NKM	Mini	50 x 28	Ivory 1/16" (1.6mm)	Non-Kit 10
10-4021-10NKP	Petite	58 x 30	Ivory 1/16" (1.6mm)	Non-Kit 10
10-4021-10NKS	Small	65 x 35	Ivory 1/16" (1.6mm)	Non-Kit 10
10-4021-10NKL	Large	69 x 43	Ivory 1/16" (1.6mm)	Non-Kit 10
10-4022-05KM	Mini	50 x 28	Blush 1/16" (1.6mm)	Kit 5
10-4022-05KP	Petite	58 x 30	Blush 1/16" (1.6mm)	Kit 5
10-4022-05KS	Small	65 x 35	Blush 1/16" (1.6mm)	Kit 5
10-4022-05KL	Large	69 x 43	Blush 1/16" (1.6mm)	Kit 5
10-4022-10NKM	Mini	50 x 28	Blush 1/16" (1.6mm)	Non-Kit 10
10-4022-10NKP	Petite	58 x 30	Blush 1/16" (1.6mm)	Non-Kit 10
10-4022-10NKS	Small	65 x 35	Blush 1/16" (1.6mm)	Non-Kit 10
10-4022-10NKL	Large	69 x 43	Blush 1/16" (1.6mm)	Non-Kit 10