


# DTR Medical<sup>®</sup> Suction Regulator



## The Design Checklist

- ✓ House style regulator with control
- ✓ Luer connector suitable for use with Verhoeven and other micro suction tubes
- ✓ Stainless steel thumb grip
- ✓ UK manufactured
- ✓ Sterile, Single-use, individually packed

### Intended Use

The house style regulator consists of a metal device with a luer connector which can be attached to various micro suction devices and tubing to regulate suction during various procedures.

### Indications for Use

It is indicated for use in ENT, General Surgery, and Orthopaedics during suctioning procedures.

It is a Single-use device, supplied sterile. This device should not be resterilised and should be disposed of once used, in accordance with accepted medical practice and applicable, local laws and regulations.

### Instructions For Use (IFU)

There is no IFU available for the Suction Regulator.

### Conformity to European Regulations

Class I Sterile device as per Rule 1, which conforms to Medical Devices Directive 93/42/EEC.

Manufactured under EN ISO:13485 in an ISO 7 Cleanroom Environment.

## DTR Medical<sup>®</sup> Suction Regulator

### Product Specification

Code	Description	Pack Size	Box Size
<b>DSRS5001</b>	Suction Regulator	1 Suction Regulator	50 Units

### Material Breakdown

Product Component	Details
<b>Suction Regulator</b>	303 Stainless Steel

